
## pH Testing Procedure:

- 1. Weigh 5.00 grams of soil directly from sampling container and place in testing vial:remove foreign material such as pebbles, roots, etc.
- 2. Add 15 ml of glass distilled water
- 3. Place magnet in vial and stir with glass rod until thoroughly mixed
- 4. Stir for 5 minutes (300 seconds) on stir plate
- 5. Place electrode in sample according to diagram below:



- 6. Stir for additional 90 seconds with electrode still in sample
- 7. Record pH during the last 15 seconds (if fluctuating, average readings during this time)
- 8. Dispose of sample and wash vial with detergent Rinse thoroughly and use glass distilled for final rinse

## Suggestions:

- 1. Weigh samples in "batch" lots (approxiamtely 20) but do not add water until just prior to pH measurement
- 2. Maintain a constant position for electrode placement; a shift in level may alter measurement
- 3. Use same water source and check pH occassionally
- 4. If samples cannot be tested for an extended period (over 1 hr) cover with Parafilm Do not use samples which have water added and have not been tested within 1 hr.
- 5. Calibrate pH machine according to manual Use two buffers Check calibration after approximately 20 samples and adjust if necessary

## Organic Matter Analysis Procedures

- 1. Weigh empty crucible and label with appropriate sample number.
- 2. Place 5.00 grams of soil in crucible-remove larger pebbles, roots, and any foreign material.
- 3. Place in drying oven at 100 degrees Celsius for 2 hours.
- 4. Reweigh sample
- 5. Place in furnace at 600 degrees Celsius for an additional 2 hours.
- 6. Reweigh sample

Suggestions:

1. Weigh out samples in batches of 16 (this is the capacity of the furnace).

 After samples are removed from each oven, allow a constant time period for cooling (ie. remove samples from oven, cool for 3 minutes and then weigh).

- 3. Use a pencil for identifying sample numbers (write on bottom of crucible).
- 4. Dispose of samples after data records have been checked Wash with detergent and dry thoroughly before next use.

Formulas:

1. Soil Moisture

$$\frac{5.00 - B}{5.00}$$
 X 100

B=weight of sample after drying at 100 degrees C

2. Organic Matter

$$\frac{B - C}{B} \times 100$$

C=weight of sample after drying at 600 degrees C